GenENDO K-Files

fr - Cathétérisme et mise en forme canalaire
en - Catherisation and root canal shaping
es - Cateterismo y conformación del conducto radicular
de - Katheterisierung und Wurzelkanalaufbereitung
it - Cateterizzazione e sagomatura del canale radicolare
nl - Katheterisatie en vormgeving van het wortelkanaal
pl - Kateteryzacja i kształtowanie kanału koronowego
cs - Katelizace a preparace kořenného kanálu
tr - Kateterizasyon ve kök kanal şekillendirme
ar - فتح مجرى سنن وتنظيف الجذرة

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9. Symboles
Non stérile
Tenir à l'abri du soleil
Conservation au sec
Assortiment le cas échéant
Quantité
Dispositif médical
Lot
Numéro de lot
K Files
Identifiant unique du dispositif
K Files
Consultez le mode d'emploi
Date de création du mode d'emploi : 18 avril 2025
Date de révision du mode d'emploi : 18 avril 2025
Instructions for use
1. Indication
Non-surgical root canal treatment:
• An irreversibly damaged or necrotic pulp
• Elective devitalization
Intended use: The medical devices are used for manual catheterism (n°008 to 015) during a non-surgical endodontic treatment. K-Files 008 and 010 are also used in the exploration, initial penetration and root canal permeabilization.
Intended user: For use by dental professionals only.
Intended patient population:
All patients, except:
• Patient which presents a contraindication below visiblement seen.
2. Contraindications
A non-surgical treatment is not indicated for:
• Teeth that cannot be made functionally or restored.
• Teeth with insufficient periodontal support.
• Teeth with poor prognosis, uncooperative patients or patients where dental treatment procedures cannot be undertaken.
• Teeth of patients with poor oral condition that cannot be improved with a reasonable process.
Precaution to be taken with patient having a risk of endocarditis heart disease.
The decision to use an endodontic instrument must be related to the expertise of the clinical case, particularly if a root canal anatomy is considered to be too complex or in case of young patient.
3. Side effects/Complications
In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, chips, stripping, zigzag, false path, perforation, etc.) could occur and lead to a risk of infectious process.
4. Characteristics and warnings
Type 1 manual instruments according to EN ISO 3630-1:2019 Materials of the operative part: Stainless steel
Non-sterile devices.
Diamètre Maximum recommended number of uses (if the file is not visibly damaged)
• Respect the good dental practice in particular by using a dental dam and gloves.
• Keep present UDI and packaging information until the last use.
• Do not use if doubtful on the storage conditions.
• In case of doubt concerning the product identification, do not use.
• Check the condition of the instrument's blade and its fit with the handle before use. If the instrument is damaged or shows signs of wear, do not use it.
• Inform the manufacturer and the national regulatory authority of any serious incident relating to the instrument.
• Be sure to respect the maximum recommended number of uses.
• Follow the reprocessing instructions for use and use and reuse.
• Medical device class according to MDR 2017/745: 1.
5. Clinical claims
There are no specific clinical claims for the devices other than the achievement of the indication for use that corresponds to paragraph 1.
6. Reprocessing instructions
• For your own safety, please wear personal protective equipment (gloves, glasses and mask).
• Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition.
• Always respect the instructions for use provided by the manufacturer of the products.
• Due to the product design and the materials used, no definite limit to the maximum number of performance processing cycles can be specified.
The service life of the medical devices is determined by their function and careful handling. Multiple use and re-sterilization cycles may lead to wear and tear of the instrument.
The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements.
The state of the art and national laws require that validated procedures be followed.
• Gloves, masks, gown, dental dam, as recommended by the cleaning agents manufacturers
• Running Tap water, deionized water or purified water according to European Pharmacopoeia
• Cleaning agent (neodisher® Septo Active and neodisher® MedZym)
• Spray soft bristles
• Container
• Ultrasonic tub or washer-disinfector
• Class B sterilization apparatus according to EN ISO 17665.
Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).

10. Limitations of representation
• For your own safety, please wear personal protective equipment (gloves, glasses and mask).
• Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition.
• Always respect the instructions for use provided by the manufacturer of the products.
• Due to the product design and the materials used, no definite limit to the maximum number of performance processing cycles can be specified.
The service life of the medical devices is determined by their function and careful handling. Multiple use and re-sterilization cycles may lead to wear and tear of the instrument.
The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements.
The state of the art and national laws require that validated procedures be followed.
• Gloves, masks, gown, dental dam, as recommended by the cleaning agents manufacturers
• Running Tap water, deionized water or purified water according to European Pharmacopoeia
• Cleaning agent (neodisher® Septo Active and neodisher® MedZym)
• Spray soft bristles
• Container
• Ultrasonic tub or washer-disinfector
• Class B sterilization apparatus according to EN ISO 17665.
Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).

11. Initial treatment
Place use products in a container or a wipe with running water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min. Rinse the products with running tap water at 20-40 °C for 5 min.
If there is time before the next step, make sure the device remains moist by placing it in a wet wipe. Do not exceed 1 hour.
Remarks:
• Do not use filing agents or hot water (40 °C), as this causes fixation of residues and can impair successful cleaning.
• Follow instructions and respect the concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
Groupe de patients previsto:
Grupo de pacientes previsto:
12. Preparation before cleaning
If there are visible impurities on a device, it is recommended to manually pre-clean by brushing under running tap water at 20-40 °C for at least 1 min with a soft brush until all impurities have been removed.
Remark: Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
13. Visual inspection
Inspect the products and discard damaged products (broken, untwisted or abnormally bent).
14. Cleaning
Insert the products into an ultrasound apparatus beaker.
Run ultrasound apparatus for 10-15 min with deionized water and 0.5-2.0% neodisher® MedZym.
Remarks: Follow instructions, observe water quality, concentrations and cleaning time specified by the manufacturer of the cleaning solution.
15. Rinsing
Rinse the products with purified water according to European Pharmacopoeia at 20-40 °C for 1 min.
Remark: It is recommended to use deionized water.
16. Drying
Dry the products with compressed air until products are visibly dry.
17. Cleaning/Rinsing/Drying
Place the instruments in the tray of the on the side-in trolley of the washer/dishwasher. Ensure that the cleaning agent and water can reach all area of the instrument's blade.
Pre-cleaning: cold tap water, 4 min
Run cleaning cycle with 0.5-1.0% neodisher® MedZym at 50°C for 10min
Final rinsing with purified water according to European Pharmacopoeia for 3 min Perform drying.
Remarks:
• Follow instructions and concentrations given by the manufacturer of the cleaning solution.
• Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer.
• Use only approved washers-disinfectors according to EN ISO 15883, maintained and validated regularly.
18. Visual inspection
Inspect the products.
Re-do steps 4-5 if the product is visibly not clean or discard any damaged product.
19. Packaging
Place the instruments in two paper-plastic pouches for steam sterilization in compliance with ISO 11607 and EN 868 standards.
Remarks:
• Seal the pouches as recommended by the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
20. Sterilization
Sterilize the products with steam vapor.
The following sterilization cycles can be used:
• Apparatus: class B
• 132 °C (269.6 °F), 4 min
• 134 °C (273.2 °F), 3 min
• 134 °C (273.2 °F), 18 min
• Minimum drying: 20 min
Control the physico-chemical indicators and cycle parameters.
We recommend a steam sterilization at 134 °C and duration of 18 min for prior inactivation according to French regulations.
Remarks:
• When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded.
• Place the pouches in the steam sterilizer according to the recommendations given by the pouch manufacturer.
• Use only a pre-vacuum air removal sterilizer that meets the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full-size sterilizer), with saturated steam.
• Use only validated sterilizer according to EN ISO 17665, maintained and validated regularly.
21. Storage
Store the products in a dry, clean and dust-free environment at the temperature specified by the product packaging with the steam sterilization manufacturer.
Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete runn should be performed.
22. Storage and transport conditions
Keep away from sunlight and keep dry.
Keep the initial packaging closed until the reprocessing step (§6).
Close immediately the initial packaging after removing an instrument.
23. Disposal
After use, instruments must be placed in a secure bin intended for cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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15. Rinçage
Rinse the products with purified water according to European Pharmacopoeia at 20-40 °C for 1 min.
Remark: It is recommended to use deionized water.
16. Drying
Dry the products with compressed air until products are visibly dry.
17. Cleaning/Rinsing/Drying
Place the instruments in the tray of the on the side-in trolley of the washer/dishwasher. Ensure that the cleaning agent and water can reach all area of the instrument's blade.
Pre-cleaning: cold tap water, 4 min
Run cleaning cycle with 0.5-1.0% neodisher® MedZym at 50°C for 10min
Final rinsing with purified water according to European Pharmacopoeia for 3 min Perform drying.
Remarks:
• Follow instructions and concentrations given by the manufacturer of the cleaning solution.
• Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer.
• Use only approved washers-disinfectors according to EN ISO 15883, maintained and validated regularly.
18. Visual inspection
Inspect the products.
Re-do steps 4-5 if the product is visibly not clean or discard any damaged product.
19. Packaging
Place the instruments in two paper-plastic pouches for steam sterilization in compliance with ISO 11607 and EN 868 standards.
Remarks:
• Seal the pouches as recommended by the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
20. Sterilization
Sterilize the products with steam vapor.
The following sterilization cycles can be used:
• Apparatus: class B
• 132 °C (269.6 °F), 4 min
• 134 °C (273.2 °F), 3 min
• 134 °C (273.2 °F), 18 min
• Minimum drying: 20 min
Control the physico-chemical indicators and cycle parameters.
We recommend a steam sterilization at 134 °C and duration of 18 min for prior inactivation according to French regulations.
Remarks:
• When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer's maximum load is not exceeded.
• Place the pouches in the steam sterilizer according to the recommendations given by the pouch manufacturer.
• Use only a pre-vacuum air removal sterilizer that meets the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full-size sterilizer), with saturated steam.
• Use only validated sterilizer according to EN ISO 17665, maintained and validated regularly.
21. Storage
Store the products in a dry, clean and dust-free environment at the temperature specified by the product packaging with the steam sterilization manufacturer.
Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete runn should be performed.
22. Storage and transport conditions
Keep away from sunlight and keep dry.
Keep the initial packaging closed until the reprocessing step (§6).
Close immediately the initial packaging after removing an instrument.
23. Disposal
After use, instruments must be placed in a secure bin intended for cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

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