# Endomethasone N/ Endomethasone Liquid

Permanent root canal sealer

## Presentation

| Presentation                                                                                            | Number of applications |
|---------------------------------------------------------------------------------------------------------|------------------------|
| 1 bottle of 14g of Endomethasone N + 1 spoon                                                            | 27*                    |
| 1 bottle of 42g of Endomethasone N + 1 spoon                                                            | 86*                    |
| 1 bottle of 10 mL Endomethasone Liquid                                                                  | 41**                   |
| 1 bottle of 14g of Endomethasone N + 1 spoon + 1 bottle of 10 mL of Endomethasone Liquid + 1 mixing pad | 27*/41**               |

\* 1 application of powder = 2 full spoons dose

\*\* 1 application = 4 to 6 drops

## Composition

### Endomethasone N:

Ivory powder:

Hydrocortisone acetate ..... 1.00g

Thymol iodide, Barium sulfate, Zinc oxide, Magnesium stearate q.s.ad ..... 100g

### Endomethasone Liquid:

Light yellow to amber clear liquid: eugenol and star anise peppermint oil.

## Intended use

Endomethasone N powder to be mixed with Endomethasone Liquid to obtain a permanent root canal sealer.

## Indications for use

Permanent root canal filling in combination with gutta-percha points in case of inflamed or necrotic pulp in the course of initial treatment or retreatment procedure.

## Target group(s)

### Target population

Permanent teeth (mature & immature)

Pregnant women and children < 30 kg: the product containing a small amount of hydrocortisone acetate, in the absence of specific studies on those populations, the use during pregnancy, during breastfeeding or in children under 30kg is possible only if the practitioner deems it necessary or in the absence of another alternative to root canal treatment.

**Intended users**

For dental healthcare professional.

**Clinical benefit**

Preservation of the tooth functionality.

Help short-term endodontic pain management.

Summary of safety and clinical performance is available, using Basic UDI-DI code 37601407000010G4, 37601407000011G6 or 37601407000011G8, in the European database on medical devices (EUDAMED): <https://ec.europa.eu/tools/eudamed> or upon request ([sscp@septodont.com](mailto:sscp@septodont.com)) until EUDAMED will be functional.

**Performances characteristics****Technical performance**

Hermetically seal of the root canal.

**Technical specifications**

| <b>Characteristics</b> | <b>Range Value<br/>According to EN ISO 6876:2012</b>         | <b>Mean Value<br/>Endomethasone N/ Endomethasone Liquid</b> |
|------------------------|--------------------------------------------------------------|-------------------------------------------------------------|
| Setting time           | According to manufacturing specifications: From 1 to 6 hours | From 1 to 6 hours<br>Average: 2.5 h $\pm$ 30min             |
| Radiopacity            | $\geq$ 3mm AI                                                | $\geq$ 6mm AI                                               |
| Working time of 5min   | $\geq$ 17mm                                                  | $\geq$ 20mm                                                 |

**Clinical performance**

Permanent root canal obturation of permanent teeth.

Help post-operative pain reduction following endodontic treatment due to the anti-inflammatory ancillary action of hydrocortisone acetate.

**Contra indication**

Hypersensitivity to one or more components.

Temporary teeth.

**Expected lifetime**

10 years.

## Instruction for use

### Mixing instructions

- 1/ Open the bottle and collect two levelled spoon of Endomethasone N powder with the spoon supplied in the box and put them on the mixing pad.
- 2/ Pour 4 to 6 drops of Endomethasone Liquid on the mixing pad next to the powder deposit.
- 3/ Add progressively the Endomethasone N powder to the Endomethasone Liquid until obtention of the desired consistency. The mix must be made extemporaneously.
- 4/ Close the bottles tightly and quickly after each use
- 5/ Immediately rinse and dry the spoon to remove any residual material.

### Placement into the root canal

**WARNING:** Perform an apexification before using the cement to avoid extrusion of filling material Do not use Endomethasone N and Endomethasone Liquid if the apexification cement isn't perfectly set. Otherwise, technical specifications won't be guaranteed and may lead to the failure of the endodontic procedure

- 1/ Perform the root canal cleaning and shaping procedure using standard endodontic procedures.
- 2/ Flush with sterile saline solution to remove any residue of irrigant and dry without totally dehydrating the root canal.
- 3/ Prepare the cement paste according to the above-mentioned mixing instructions.
- 4/ Apply a coating of cement paste onto the canal walls using either a paper point or a gutta percha point.
- 5/ Complete the obturation by inserting one or several gutta-percha point(s).

### Retreatment procedure

Remove the existing sealer using conventional removal techniques. Then follow instructions given in the sections "Mixing instructions" and "Placement into the root canal".

## Warnings and precautions for use

### Precautions for use

- As Endomethasone N contains hydrocortisone acetate, it is at the discretion of the practitioner to choose the product according to the patient's clinical condition and current endodontic practices.
- After first opening:
  - use Endomethasone N 14 g within 6 months
  - use Endomethasone N 42 g within 9 months
  - use Endomethasone Liquid within 9 months
- The device may cause eyes and skin irritation. Do not manipulate the device without appropriate personal protective equipment.
- Use a rubber dam to ensure that the cement is not in contact with the mucosa or swallowed.
- Technical specifications are not guaranteed and may lead to the failure of the endodontic procedure if the following is not respected:
  - . Use the spoon provided with Endomethasone N.
  - . Do not use Endomethasone N with another liquid than Endomethasone Liquid and vice versa.
  - . Use the recommended ratio of liquid/powder.

## Warnings

- Do not use the paste on several patients to avoid cross contamination.
- The working length in the root canal should not exceed 2 mm from the radiological apex to avoid apical extrusion.
- Do not use in combination with polymer based composite materials, since eugenol-based preparations inhibit their polymerization.
- Do not use the cement without gutta-percha points as this will jeopardize endodontic and retreatment procedures.

## Interactions with other drugs

The amount of hydrocortisone acetate used in root canal treatment is not likely to cause medically significant drug interactions.

## Side-effects

As part of the endodontic procedure, residual pain and inflammation could occurred.  
Possible hypersensitivity to one of the components of Endomethasone N and/or Endomethasone Liquid.  
Please report any serious incident that occurred in relation to the device to Septodont ([vigilance@septodont.com](mailto:vigilance@septodont.com)) and to the local competent authority for medical devices.

## Disposal

Endomethasone N is toxic to aquatic life with long lasting effects.  
After use put the cement and bottles in an appropriate chemical waste bin.

## Storage


Do not store above 30°C; protect from moisture.

## Access to the instructions for use

Please check regularly to see if the instructions for use have not been updated.  
The latest electronic instruction for use as well as historical versions can be downloaded and printed directly from the website <https://www.septodontcorp.com/endomethasone-n-liquid-eifu/>.

A paper copy can be requested through the website <https://www.septodontcorp.com/endomethasone-n-liquid-eifu/>.

## Symbols

Period-after-opening

Authorised Representative in CH

UK responsible person

UK marking

## EU-IFU must contain the following information:

- Date of CE mark:  
(CE): YYYY-MM
- Date of latest revision of the IFU:  
YYYY-MM (to be determined when the artwork will be done)
- IFU reference:  
S05 xxx xxx xx xx
- Manufacturer & manufacturer address:

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- UKCA marking

- UKREP:

Callisto Regulatory Consulting Limited, Iron Farm, 7  
Grimes Gate, Diseworth, Derby, DE74 2QD, United  
Kingdom

- CH REP:

SMQS  
by Arnold Bott AG  
Industriestrasse 59  
8152 Glattbrugg

# Endomethasone N/ Endomethasone Liquid

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